


### Lens Geometry


Base Curve	Diameter (mm)	Decenter (mm)	Effective Diameter (mm)	True Curve (1.53)	Convex Radius (mm)	Back Curve (1.53)	Center Thickness (mm)	Edge Thickness (mm)
2.50	75	2.5	80	2.62	202.45	2.44	12.13	11.28
4.00	75	2.5	80	4.16	127.27	5.57	8.28	9.98
5.00	75	2.5	80	5.19	102.09	5.42	9.27	9.40

### Material

- Composite 59 Polycarbonate
- $n_d = 1.586$
- Abbe 30
- Specific Gravity 1.20 <sup>g</sup>/<sub>ccm</sub>
- 100% UV-A and UV-B Absorption



**Right Lens**



**Left Lens**

Permanent Engravings	Nasal Engraving	Temple Engraving
<b>atLast! Composite 59</b>	○ !59	○ 250
<b>atLast! Composite 59 Transitions® Lenses (Gray &amp; Brown)</b> <i>(Launching July, 2009)</i>	!59T	250

**Minus Cylinder**

	0.00	-0.25	-0.50	-0.75	-1.00	-1.25	-1.50	-1.75	-2.00	-2.25	-2.50	-2.75	-3.00	-3.25	-3.50	-3.75	-4.00
3.50																	
3.25																	
3.00																	
2.75																	
2.50																	
2.25																	
2.00																	
1.75																	
1.50																	
1.25																	
1.00																	
0.75																	
0.50																	
0.25																	
0.00																	
-0.25																	
-0.50																	
-0.75																	
-1.00																	
-1.25																	
-1.50																	
-1.75																	
-2.00																	
-2.25																	
-2.50																	
-2.75																	
-3.00																	
-3.25																	
-3.50																	
-3.75																	
-4.00																	
-4.25																	
-4.50																	
-4.75																	
-5.00																	
-5.25																	
-5.50																	
-5.75																	
-6.00																	

**Prescription Range**

**(Composite 59 Clear and Composite 59 Transitions Gray or Brown – Transitions<sup>®</sup> Launching in mid July)**

- Sphere: +3.50D to -6.00D
- Cylinder: up to -4.00D
- Add power: 1.00D to 2.50D

**www.atlastlens.com • 888-749-3567**

**U.S. Patent No. 7,475,985. Other US and Foreign Patents Pending.**